



DWG NO: AHE370

FOR REFERENCE ONLY

NOTES:

- 1 MATERIAL: PVDF
USP CLASS VI
RoHS COMPLIANT
FDA 21 CFR 177.2510
ANIMAL FREE
COLOR: NATURAL
- 2 ALL DIMENSIONS ARE TYPICAL



SCALE 3:1

 CPC™ <i>Smart fluid handling to take you forward, faster</i> www.cpcworldwide.com	
PROJECT: 9975	3/32 HB ELBOW UNION (A-SERIES) (CUSTOMER DRAWING)
DWG BY: D.CHEN	DATE: Jan-16-18 CHKD BY:
MAT'L: SEE NOTE 1	DWG NO: AHE370