



DWG NO: ALF41

FOR REFERENCE ONLY

NOTES:

☐ MATERIAL: POLYPROPYLENE
FDA 21 CFR 177.1520
RoHS COMPLIANT
ANIMAL FREE
USP CLASS VI
COLOR: NATURAL



 Smart fluid handling to take you forward, faster www.cpcworldwide.com	
PROJECT: I-1005	1/8 HB FEMALE LUER (CUSTOMER DRAWING)
DWG BY: D.CHEN	DATE: 06-29-17
MAT'L: SEE NOTE ☐	CHKD BY:
	DWG NO: ALF41