

DWG NO: HER43


FOR REFERENCE ONLY

NOTES:

☐ MATERIAL: POLYPROPYLENE  
FDA 21 CFR 177.1500  
RoHS COMPLIANT  
ANIMAL FREE  
USP CLASS VI  
COLOR: NATURAL



SCALE 3:1



 **CPC™** *Smart fluid handling to take you forward, faster*  
[www.cpcworldwide.com](http://www.cpcworldwide.com)

PROJECT:	1/8 HB X 3/32 HB REDUCING ELBOW (CUSTOMER DRAWING)		
7362			
DWG BY: GJH	DATE: Jan-16-18	CHKD BY:	
MAT'L: SEE NOTE ☐	DWG NO: HER43		