

DWG NO: HSR84-CUST


CONFIDENTIAL

**THIS PRINT CONTAINS PROPRIETARY INFORMATION,
NO REPRODUCTION OR DISCLOSURE OF THIS DOCUMENT
SHALL BE MADE WITHOUT WRITTEN AUTHORIZATION FROM
COLDER PRODUCTS COMPANY**

NOTES:

MATERIAL: POLYPROPYLENE
FDA 21 CFR 177.1500
ROHS COMPLIANT
ANIMAL FREE
USP CLASS VI

COLOR : NATURAL



PARTING LINE FREE HOSE BARB

		Colder Products Company®	
		www.colder.com 651 645-0091 1 800 444-2474	
PROJECT	1/4 X 1/8 REDUCER		
7362	(CUSTOMER DRAWING)		
DWG BY: VS	DATE: 07/20/09	CHKD BY:	
MATL/FIN: SEE NOTE <input type="checkbox"/>		DWG NO: HSR84-CUST	