



DWG NO: LF3131-CUST

NOTES:

1 MATERIAL: NYLON
FDA 21 CFR 177.1500
RoHS COMPLIANT
COLOR: BLACK



| | | | |
|----------------------|---|---------------------|------------------|
| CPC | | 1 800 444 2474 | cpcworldwide.com |
| PROJECT 7362 | 3/32 HB FEMALE LUER (CUSTOMER DRAWING) | | |
| DWG BY: KRN | DATE: 07/01/09 | CHKD BY: | |
| MATL/FIN: SEE NOTE 1 | | DWG NO: LF3131-CUST | |