


DWG NO: LFT30

FOR REFERENCE ONLY

NOTES:

□ MATERIAL: NYLON
FDA 21 CFR 177.1500
RoHS COMPLIANT
COLOR: WHITE



SCALE 2:1

 **CPC**™ *Smart fluid handling to take you forward, faster*
www.cpcworldwide.com

PROJECT:	FEMALE TEE LUER (CUSTOMER DRAWING)	
8530		
DWG BY: SOMA	DATE: Jan-16-18	CHKD BY:
MAT'L: SEE NOTE □	DWG NO: LFT30	