


DWG NO: PMLF4131-CUST

NOTES:

- 1 MATERIAL: NYLON
FDA 21 CFR 177.1500
RoHS COMPLIANT
COLOR: BLACK
- 2 PANEL NUT (PMNUTSST) AND COLOR-CODE RING (PMRLXX)
ARE AVAILABLE SEPERATELY



1 800 444 2474

cpcworldwide.com

PROJECT
7362

1/8 HB PANEL MOUNT FEMALE LUER
(CUSTOMER DRAWING)

DWG BY: AMK

DATE: 07/06/09

CHKD BY:

MATL/FIN: SEE NOTE 1

DWG NO: PMLF4131-CUST