


DWG NO: SLM6130-CUST

NOTES:

- 1 MATERIAL: NYLON
FDA 21 CFR 177.1500
RoHS COMPLIANT
COLOR: WHITE
- 2 CAN BE USED WITH LMSRXX (ROTATING LUER LOCK RING) OR LMSLXX (STATIONARY LUER LOCK RING)



1 800 444 2474

cpcworldwide.com

PROJECT
7362

3/16 HB MALE SLIP LUER
(CUSTOMER DRAWING)

DWG BY: PRA

DATE: 06-30-09

CHKD BY:

MATL/FIN: SEE NOTE 1

DWG NO: SLM6130-CUST