


DWG NO: SLMU

FOR REFERENCE ONLY

NOTES:

- MATERIAL: POLYPROPYLENE
 FDA 21 CFR 177.1520
 RoHS COMPLIANT
 ANIMAL FREE
 USP CLASS VI
 COLOR: NATURAL
- ALL DIMENSIONS ARE TYPICAL



CPC™ *Smart fluid handling to take you forward, faster*
www.cpcworldwide.com

PROJECT:	FITTING, MALE SLIP X MALE SLIP UNION NAT (CUSTOMER DRAWING)	
9553		
DWG BY: R.SI	DATE: Jan-16-18	CHKD BY:
MAT'L: SEE NOTE <input type="checkbox"/>	DWG NO: SLMU	